

RETAINER BRITE™

Professional Retainer & Appliance Cleaning System

Downloaded from: <https://theretainerbrite.com/>

INSTRUCTIONS FOR USE

What is Retainer Brite?

Retainer Brite is a professional-grade cleaning system specifically formulated to clean and freshen retainers, aligners, mouthguards, and other dental appliances. The powerful oxygenating formula removes plaque, bacteria, and odors while being gentle on appliance materials.

How to Use Retainer Brite

Daily Cleaning Method:

1. **Fill container** with warm (not hot) water - approximately 3-4 ounces
2. **Drop one tablet** into the water and allow it to dissolve completely
3. **Place your appliance** into the solution immediately after tablet dissolves
4. **Soak for 15 minutes** for daily cleaning
5. **Remove appliance** and rinse thoroughly with clean water
6. **Brush gently** if needed with a soft-bristled toothbrush
7. **Rinse again** before placing back in your mouth

Deep Cleaning Method:

- For heavily soiled appliances or weekly deep cleaning
- Follow steps 1-3 above
- **Extend soaking time to 30 minutes**
- Complete steps 5-7 as normal

Important Safety Information

Do's:

- Use warm water (never hot or boiling)
- Allow tablet to completely dissolve before adding appliance
- Rinse appliance thoroughly after soaking
- Store tablets in a cool, dry place

- Keep container tightly closed when not in use

Don'ts:

- Do not use hot water - may damage appliance materials
- Do not exceed recommended soaking time
- Do not swallow the cleaning solution
- Do not use on metal appliances with electrical components
- Do not store dissolved solution for later use

Frequency of Use

- **Daily users:** Clean every day after removal
- **Part-time users:** Clean after each use
- **Storage cleaning:** Clean before storing appliance for extended periods

Troubleshooting

White residue on appliance:

- Rinse more thoroughly after cleaning
- Use cooler water temperature

Strong taste after cleaning:

- Increase rinse time with clean water
- Brush gently with water only

Appliance discoloration:

- Normal for heavily stained appliances
 - May require multiple treatments
-

INGREDIENTS

Active Ingredients:

- **Sodium Bicarbonate** - Gentle abrasive and pH buffer
- **Citric Acid** - Natural cleaning agent and freshener
- **Sodium Carbonate Peroxide** - Oxygenating bleaching agent
- **Potassium Monopersulfate** - Antimicrobial oxidizing agent

Inactive Ingredients:

- **Sodium Lauryl Sulfoacetate** - Gentle foaming agent
 - **Polyethylene Glycol** - Tablet binding agent
 - **Sodium Benzoate** - Preservative
 - **Natural Mint Flavoring** - Freshening agent
 - **Blue #1 Lake** - Coloring agent
-

STORAGE AND SAFETY

Storage Instructions:

- Store in original container at room temperature
- Keep away from moisture and direct sunlight
- Keep out of reach of children and pets
- Do not transfer to unmarked containers

Safety Precautions:

- **For external use only** - not for ingestion
- **Keep away from eyes** - if contact occurs, flush with water
- **Wash hands** after handling tablets
- **Discontinue use** if irritation occurs
- **Consult dentist** if you have concerns about appliance cleaning

Emergency Information:

- **If swallowed:** Rinse mouth with water and contact poison control
 - **If eye contact:** Flush immediately with clean water for 15 minutes
 - **If skin irritation:** Discontinue use and wash affected area
-

FREQUENTLY ASKED QUESTIONS

Q: Can I use Retainer Brite on all types of retainers? A: Yes, Retainer Brite is safe for most retainer materials including plastic, acrylic, and wire retainers. Not recommended for electronic or battery-powered devices.

Q: How often should I replace my cleaning solution? A: Always use fresh solution. Do not reuse dissolved cleaning solution.

Q: Can I use this product for dentures? A: While safe for most dental appliances, consult your dentist for denture-specific recommendations.

Q: What if my appliance has metal parts? A: Retainer Brite is safe for most metal components found in retainers and braces.

CONTACT INFORMATION

For questions about this product or to report adverse effects:

Website: <https://theretainerbrite.com/>

Customer Service: Available through website contact form

Product Support: Online FAQ and troubleshooting guides available

This document was downloaded from: <https://theretainerbrite.com/>

Retainer Brite™ - Professional Appliance Care

Keep your smile bright and your appliances clean

Revision Date: May 2025

Document Version: 1.0